

RTqPCR: Functions for analysis and normalisation of pre-processing RT-qPCR data

Navneet Phogat and Matthias Kohl
Institute of Precision Medicine
Hochschule Furtwangen University, Germany

May 16, 2017

Contents

1	Introduction	1
2	Reading the raw fluorescent data	2
3	Reading the sample information data	4
4	Computing the Cq (Cycle Threshold) and amplification efficiencies	6
5	Replacing the values and removing the non-detects of Cq values and amplification efficiencies	11
5.1	ReplaceAboveCutOff	11
5.2	ReplaceValue	12
5.3	ReplaceNAs	13
5.4	NonDetects	14
6	Combining the replicates, based on Cq values and amplification efficiencies	15
7	Compute delta Cq, delta delta Cq values and fold concentration of the Target samples	21
7.1	Compute delta Cq values	21
7.2	Compute delta delta Cq values and fold concentration	22
8	Compute relative expression ratio	23

9	Auxiliary functions	23
9.1	<code>read.Mx3005P</code>	23
9.2	<code>Read.LC480SampleInfo</code>	25
9.3	<code>Read.Mx3005PSampleInfo</code>	25
9.4	<code>DeltaDeltaCq</code>	26
9.5	<code>NRQeffs</code>	27
9.6	<code>RTqPCR.dataframe</code>	27
10	Graphical User Interface (GUI)	27

1 Introduction

The package "RTqPCR" provides methods for the normalization and analysis of pre-processing real-time quantitative RTqPCR data. In this vignette we describe and demonstrate the available functions. The package "RTqPCR" is designed on the basis of experimental data of LC480 light cycler and Mx3005P RT-qPCR. If the user wants to implement the package for any other RTqPCR, then, user may change the data into one of the two formats of LC480 light cycler and Mx3005P RTqPCR. The package "RTqPCR" is divided into six parts: First, we show how the user may read the fluorescence and sample information raw data files and compute Cq values and amplification efficiencies. Second, we demonstrate, how the user may deal with Cq values and efficiencies above user-chosen threshold values, replace the particular values and deal with non-detects in Cq values and efficiencies. Third, how the user may combine technical replicates and compute standard deviation of Cq values and amplification efficiencies within replicates. Fourth, how the user may compute delta Cq values and standard deviation of delta Cq, compute delta delta Cq values, standard deviation of delta delta Cq values, fold concentration (log values) based on delta delta Cq values and visualisation of the results of fold concentration in the form of bar plots too. Fifth, how the user may compute the relative expression ratio based on the three different methods - Method by Roche, Method by Pfaffl et al., ddCq (delta delta Cq) method and may visualise the results in the form of bar plots too. Sixth part is about auxiliary functions.

User may load the "RTqPCR" package.

```
> library(RTqPCR)
```

2 Reading the raw fluorescent data

The `read.RTqPCR` function is based on `[ReadqPCR]read.LC480` function from the "Read-qPCR" package from Bioconductor software in R software and `read.Mx3005P` function. Details of `read.LC480` function is present in "ReadqPCR" package. The detail of `read.Mx3005P` is present in section Auxiliary Function. The `read.RTqPCR` function reads in the .txt Tab delimited file of raw fluorescent data of both LC480 light cycler and Mx3005P RTqPCR

and populate an object of class `RTqPCRBatch` with new slots. Following is the example to read in the raw fluorescent data of LC480 light cycler through `read.RTqPCR` function:

```
> #library(ReadqPCR) # load the ReadqPCR library
> path <- system.file("exData", package = "RTqPCR")
> LC480.example <- file.path(path, "LC480_Example.txt")
> rtData.LC480 <- read.RTqPCR(LC480.example, PCRtype = "LC480")
> rtData.LC480 ## to express the overview of fluorescent data
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 45 features, 96 samples
  element names: exprs
protocolData: none
phenoData
  sampleNames: A1 A2 ... H12 (96 total)
  varLabels: Sample position Sample name Program number Segment number
  varMetadata: labelDescription
featureData
  featureNames: 1 2 ... 45 (45 total)
  fvarLabels: Cycle number Acquisition time Acquisition temperature
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> head(exprs(rtData.LC480[,1:5])) ## to express the first five fluorescent data
```

	A1	A2	A3	A4	A5
1	15.58	16.93	27.32	16.60	16.88
2	15.53	16.96	27.69	16.68	16.91
3	15.49	16.93	27.56	16.67	16.95
4	15.47	16.97	27.62	16.63	16.87
5	15.47	16.88	27.56	16.65	16.92
6	15.41	16.90	27.55	16.65	16.87

```
> head(pData(rtData.LC480)) ## to express the pheno data
```

	Sample position	Sample name	Program number	Segment number
A1	A1	Sample_1	2	3
A2	A2	Sample_2	2	3
A3	A3	Sample_2	2	3
A4	A4	Sample_3	2	3
A5	A5	Sample_4	2	3
A6	A6	Sample_4	2	3

```
> head(fData(rtData.LC480)) ## to express the feature data
```

	Cycle number	Acquisition time	Acquisition temperature
1	1	691600	71.79
2	2	770133	71.64
3	3	848716	71.64
4	4	927233	71.67
5	5	1005816	71.67
6	6	1084316	71.67

The example to read in the raw fluorescent data of Mx3005P RTqPCR through `read.RTqPCR` function is as follows:

```
> Mx3005P.example <- file.path(path, "Mx3005P_Example.txt")
> rtData.Mx <- read.RTqPCR(Mx3005P.example, PCRtype = "Mx3005P")
> rtData.Mx ## to express the overview of fluorescent data
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 45 features, 34 samples
  element names: exprs
protocolData: none
phenoData
  sampleNames: 14 15 ... 78 (34 total)
  varLabels: Well Ramp/Plateau ... mspSegment (5 total)
  varMetadata: labelDescription
featureData
  featureNames: 1 2 ... 45 (45 total)
  fvarLabels: Cycle# Temperature
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> head(exprs(rtData.Mx[,1:5])) ## to express the first five fluorescent data
```

	14	15	16	17	18
1	4827	4910	4900	4874	5127
2	4708	4801	4780	4766	5003
3	4662	4762	4721	4711	4934
4	4633	4757	4705	4688	4933
5	4634	4748	4708	4697	4928
6	4615	4749	4705	4693	4939

```
> head(pData(rtData.Mx)) ## to express the pheno data
```

	Well	Ramp/Plateau	Ramp/Plateau#	Dye	mspSegment
14	14	P	2	Campy	2
15	15	P	2	Campy	2
16	16	P	2	Campy	2
17	17	P	2	Campy	2
18	18	P	2	Campy	2
20	20	P	2	Campy	2

```
> head(fData(rtData.Mx)) ## to express the feature data
```

	Cycle#	Temperature
1	1	59.6
2	2	59.6
3	3	59.6
4	4	59.6
5	5	59.6
6	6	59.6

For other RTqPCR machines, the user needs to save the raw fluorescent data into the provided .txt Tab delimited format files of LC480 or Mx3005P.

3 Reading the sample information data

`read.RTqPCRSampleInfo` function is based on `Read.LC480SampleInfo` and `Read.Mx3005PSampleInfo` functions. The `read.RTqPCRSampleInfo` function populates sample information data into the resulting `AnnotatedDataFrame` format. The example to read in the .txt Tab delimited files of sample information data of LC480 light cycler through `read.RTqPCRSampleInfo` function, is as follows:

```
> SampleInfoLC480 <- file.path(path, "LC480_example_SampleInfo.txt")
> samInfoLC480 <- read.RTqPCRSampleInfo(SampleInfoLC480, PCRtype = "LC480")
> samInfoLC480 ## to express the overview of the sample information data
```

```
An object of class 'AnnotatedDataFrame'
 rowNames: 1 2 ... 96 (96 total)
 varLabels: Sample position Sample name ... Combined sample and target
            type (9 total)
 varMetadata: labelDescription
```

```
> head(pData(samInfoLC480))
```

	Sample position	Sample name	Replicate of	Filter combination	Target name
1	A1	Sample_1	A1	465-510	negativ Kontrolle
2	A2	Sample_2	A2	465-510	goi
3	A3	Sample_2	A3	465-510	hk
4	A4	Sample_3	A1	465-510	negativ Kontrolle
5	A5	Sample_4	A2	465-510	goi
6	A6	Sample_4	A3	465-510	hk

	Sample Pref color	Concentration	Efficiency	Combined sample and target type
1	\$00FF8000	NA	2	Target Negative
2	clRed	NA	2	Target Unknown
3	\$0030D700	NA	2	Ref Unknown
4	clFuchsia	NA	2	Target Negative
5	clGray	NA	2	Target Unknown
6	\$0012D7FA	NA	2	Ref Unknown

To read in the .txt Tab delimited file of sample information data of Mx3005P, designed on the basis of given sample format in exData of "RTqPCR" package, the example is illustrated below:

```
> SampleInfoMx <- file.path(path, "Mx3005P_example_SampleInfo.txt")
> samInfoMx <- read.RTqPCRSampleInfo(SampleInfoMx, PCRtype = "Mx3005P")
> samInfoMx ## to express the overview of the sample information data
```

An object of class 'AnnotatedDataFrame'

rowNames: 1 2 ... 34 (34 total)

varLabels: Well Replicate of Target name Combined sample and target
type

varMetadata: labelDescription

```
> head(pData(samInfoMx))
```

	Well	Replicate of	Target name	Combined sample and target type
1	14	A1	goi	Target Unknown
2	15	A2	goi	Target Unknown
3	16	A3	hk	Ref Unknown
4	17	A1	goi	Target Unknown
5	18	A2	goi	Target Unknown
6	20	A3	hk	Ref Unknown

4 Computing the Cq (Cycle Threshold) and amplification efficiencies

The `CqEfts` function compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies for LC480 light cycler and Mx3005P RTqPCR from raw fluorescent data. The `CqEfts` function works on object of class `RTqPCRBatch` and populate the final results into an object of class `RTqPCRBatch` with new slots. `CqEfts` function is based on the function `[qpcR]pcrbatch` from "qpcR". The default method of computation is the sigmoidal model. To compute the Cq values and amplification efficiencies, first of all we need to merge the fluorescent and sample information data. Following is the example to merge the fluorescent and sample information data of LC480 light cycler:

```
> ## To merge the fluorescent and sample information data
> merge.LC480<-merge(rtData.LC480,samInfoLC480)
> ## to express the overview of merge data
> merge.LC480
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 45 features, 96 samples
  element names: exprs
protocolData: none
phenoData
  sampleNames: 1 2 ... 96 (96 total)
  varLabels: Sample position Sample name ... Combined sample and target
             type (11 total)
  varMetadata: labelDescription
featureData
  featureNames: 1 2 ... 45 (45 total)
  fvarLabels: Cycle number Acquisition time Acquisition temperature
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

Following is the example to compute the CqValues, amplification efficiencies and standard errors of amplification efficiencies of LC480 light cycler data by default sigmoidal method:

```
> ## To compute the CqValues and amplification efficiencies by sigmoidal model
> res.LC <- CqEfts(merge.LC480, PCRtype = "LC480", baseline = "none")

> #to see the overview of data
> res.LC
```

```

RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 96 features, 1 samples
  element names: effs, exprs, se.effs
protocolData: none
phenoData: none
featureData
  featureNames: 1 *2* ... 96 (96 total)
  fvarLabels: Sample position Sample name ... sig.model (35 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

```

```

> #to express the first five CqValues
> exprs(res.LC)[1:5]

```

```

[1] 45.00    NA 24.87 29.37    NA

```

```

> #to express the first five amplification efficiencies
> effs(res.LC)[1:5]

```

```

[1] 1.0024280          NA 0.6576489 1.0003550          NA

```

```

> #to express the first five standard errors of amplification efficiencies
> se.effs(res.LC)[1:5]

```

```

[1] 4.560011e-02          NA 6.429260e+03 3.701013e-02          NA

```

To express all the Cq values, amplification efficiencies and standard errors in amplification efficiencies, user can implement the following code:

```

> ## To express all Cq values
> exprs(res.LC)
> ## To express all amplification efficiencies
> effs(res.LC)
> ## To express standard errors in amplification efficiencies
> se.effs(res.LC)

```

To compute the Cq values, amplification efficiencies and standard error of amplification efficiencies of LC480 light cycler data by fit exponential method:

```

> res.LC1<- CqEffs(merge.LC480, PCRtype = "LC480", Effmethod = "expfit",baseline = "none")
> res.LC1
> exprs(res.LC1)[1:5]

```



```

> effs(res.LC1)[1:5]
> se.effs(res.LC1)[1:5]
> exprs(res.LC1)
> effs(res.LC1)
> se.effs(res.LC1)

```

To compute the Cq values, amplification efficiencies and standard error of amplification efficiencies of LC480 light cycler data by window of linearity method:

```

> res.LC2<- CqEfs(merge.LC480, PCRtype = "LC480", Effmethod = "sliwin", baseline = "none")
> res.LC2
> exprs(res.LC2)[1:5]
> effs(res.LC2)[1:5]
> se.effs(res.LC2)[1:5]
> exprs(res.LC2)
> effs(res.LC2)
> se.effs(res.LC2)

```

To compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies of LC480 light cycler data by linear regression of efficiency method:

```

> res.LC3<- CqEfs(merge.LC480, PCRtype = "LC480", Effmethod = "LRE", baseline = "none")
> res.LC3
> exprs(res.LC3)[1:5]
> effs(res.LC3)[1:5]
> se.effs(res.LC3)[1:5]
> exprs(res.LC3)
> effs(res.LC3)
> se.effs(res.LC3)

```

To compute the Cq values and amplification efficiencies for Mx3005P RTqPCR, user needs to implement the "Mx3005" in argument `PCRtype` of `CqEfs` function. Remaining arguments will remain same as implemented above for computation of Cq values and amplification efficiencies for LC480 light cycler data. For example, like LC480 light cycler, here we first merge the fluorescent and sample information data of Mx3005P RTqPCR and then, compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies for Mx3005 RTqPCR. Following is the example to merge the data:

```

> ## To merge the fluorescent and sample information data
> merge.Mx<-merge(rtData.Mx,samInfoMx)

```

Following is the example to compute Cq values and amplification efficiencies for Mx3005P RTqPCR by default sigmoidal model:

```

> ## To compute the CqValues and amplification efficiencies
> ##by default method (sigmoidal model)
> res.Mx <- CqEfs(merge.Mx, PCRtype = "Mx3005P", baseline = "none")

> #To express the overview of data
> res.Mx

```

```

RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 34 features, 1 samples
  element names: effs, exprs, se.effs
protocolData: none
phenoData: none
featureData
  featureNames: 1 2 ... 34 (34 total)
  fvarLabels: Well Ramp/Plateau ... sig.model (32 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

```

```

> #to express the first five CqValues
> exprs(res.Mx)[1:5]

```

```

[1] 35.39 30.93 26.17 23.10 19.79

```

```

> #to express the first five amplification efficiencies
> effs(res.Mx)[1:5]

```

```

[1] 1.097644 1.093488 1.100702 1.108600 1.104106

```

```

> #to express the first five standard errors of efficiencies
> se.effs(res.Mx)[1:5]

```

```

[1] 34.12339 30.02550 44.42449 39.89076 37.30819

```

To visualise all the CqValues and amplification efficiencies and standard errors of amplification efficiencies, user can implement following code:

```

> ## To express all Cq values
> exprs(res.Mx)
> ## To express all amplification efficiencies
> effs(res.Mx)
> ## To express all standard errors in amplification efficiencies
> se.effs(res.Mx)

```

Following is the example to compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies for Mx3005P RTqPCR by fit exponential method:

```
> ## To compute the CqValues and amplification efficiencies by fit exponential method
> res.Mx1 <- CqEfs(merge.Mx, PCRtype = "Mx3005P", Effmethod = "expfit", baseline = "none")
> res.Mx1
> exprs(res.Mx1)[1:5]
> effs(res.Mx1)[1:5]
> se.effs(res.Mx1)[1:5]
> exprs(res.Mx1)
> effs(res.Mx1)
> se.effs(res.Mx1)
```

Following is the example to compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies for Mx3005P RTqPCR by window of linearity method:

```
> # To compute the CqValues and amplification efficiencies by window of linearity (sliwin)
> res.Mx2 <- CqEfs(merge.Mx, PCRtype = "Mx3005P", Effmethod = "sliwin", baseline = "none")
> res.Mx2
> exprs(res.Mx2)[1:5]
> effs(res.Mx2)[1:5]
> se.effs(res.Mx2)[1:5]
> exprs(res.Mx2)
> effs(res.Mx2)
> se.effs(res.Mx2)
```

Following is the example to compute the Cq values, amplification efficiencies and standard errors of amplification efficiencies for Mx3005P RTqPCR by linear regression of efficiency method:

```
> res.Mx3 <- CqValues(merge.Mx, PCRtype = "Mx3005P", Effmethod = "LRE", baseline = "none")
> res.Mx3
> exprs(res.Mx3)[1:5]
> effs(res.Mx3)[1:5]
> se.effs(res.Mx3)[1:5]
> exprs(res.Mx3)
> effs(res.Mx3)
> se.effs(res.Mx3)
```

5 Replacing the values and removing the non-detects of Cq values and amplification efficiencies

This part is further divided into four different subsections, based on the functions. The subsections are described in detail as follows:

5.1 ReplaceAboveCutOff

The function `ReplaceAboveCutOff` is specifically designed to replace the Cq values and amplification efficiencies above a cut off level by NA. The cut off level for Cq values and efficiencies is defined by the user based on the values and his experience with analysis of RT-qPCR data. Even then, it is recommended that the user should set the Cq cut off as a non-negative value below 40 and amplification efficiency cut off as non negative value below 2. The user can define the cutoff value for Cq and amplification efficiency based on his experimental data. The `ReplaceAboveCutOff` function uses the object of class `RTqPCRBatch` and returns the resulting data as an object of class `RTqPCRBatch`. The output of `CqEfs`, `ReplaceNAs`, `ReplaceValue`, `NonDetects` and `ReplaceAboveCutOff` can be implemented as an input for the `ReplaceAboveCutOff` function. Following is the example to implement the `ReplaceAboveCutOff` function:

```
> Cqeefs.cutoff <- ReplaceAboveCutOff(res.LC, NewVal = NA, Cqcutoff = 37, effscutoff = 2)
> ## to visualise the overview of the resulting data
> Cqeefs.cutoff

RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 96 features, 1 samples
  element names: efs, exprs, se.efs
protocolData: none
phenoData: none
featureData
  featureNames: 1 *2* ... 96 (96 total)
  fvarLabels: Sample position Sample name ... sig.model (35 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

> ## to express the first five Cq values
> exprs(Cqeefs.cutoff)[1:5]

[1]    NA    NA 24.87 29.37    NA

> ## to express the first five amplification efficiencies
> efs(Cqeefs.cutoff)[1:5]
```

```
[1] 1.0024280      NA 0.6576489 1.0003550      NA
```

5.2 ReplaceValue

The `ReplaceValue` function provides an opportunity to replace any value of Cq and amplification efficiency by a desired value, including NA too. This function is very useful, when because of the experimental error, manual error, instrumental error or software error the Cq and amplification efficiency for a negative sample is produced. In this case, the user can replace the values of the specific negative sample by non detects (NA). This function is also of greater use, if the user wants to replace any specific non detects by a value, may be by the value of another perfectly expressed technical replicate. The `ReplaceValue` function works on an object of class `RTqPCRBatch` and returns the resulting data too in the format of class `RTqPCRBatch`. The input for the `ReplaceValue` function can be output of any of the functions, including `CqEffs`, `ReplaceNAs`, `NonDetects`, `ReplaceAboveCutOff`, `ReplaceValue`. Following is the example of implementation of `ReplaceValue` function:

```
> Cqeffs.replace <- ReplaceValue(res.LC, NewCq = 36, Neweffs = 1, Cqrow = 24, effsrow = 24)
> ## to visualise the overview of the resulting data
> Cqeffs.replace
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 96 features, 1 samples
  element names: effs, exprs, se.effs
protocolData: none
phenoData: none
featureData
  featureNames: 1 *2* ... 96 (96 total)
  fvarLabels: Sample position Sample name ... sig.model (35 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> ## to visualise the replaced value of Cq
> exprs(Cqeffs.replace)[24]
```

```
[1] 36
```

```
> ## to visualise the replaced value of amplification efficiency
> effs(Cqeffs.replace)[24]
```

```
[1] 1
```

5.3 ReplaceNAs

The `ReplaceNAs` function is designed for replacing the non detects by a defined value. The defined value depends on the user. The `ReplaceNAs` function works on the object of class `RTqPCRBatch` and returns the resulting data too in the same format of class `RTqPCRBatch`. The output of any of the functions, including `CqEffs`, `ReplaceAboveCutOff`, `ReplaceValue`, `NonDetects` can be implemented as an input for the `ReplaceNAs` function. Following is an example to execute `ReplaceNAs` function:

```
> Cqeffs.NA <- ReplaceNAs(res.LC, NewCqNA = 30, NeweffsNA = 1)
> ## to visualise the overview of the resulting data
> Cqeffs.NA

RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 96 features, 1 samples
  element names: effs, exprs, se.effs
protocolData: none
phenoData: none
featureData
  featureNames: 1 *2* ... 96 (96 total)
  fvarLabels: Sample position Sample name ... sig.model (35 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

> ## to express first five Cq values
> exprs(Cqeffs.NA)[1:5]

[1] 45.00 30.00 24.87 29.37 30.00

> ## to express first five amplification efficiencies
> effs(Cqeffs.NA)[1:5]

[1] 1.0024280 1.0000000 0.6576489 1.0003550 1.0000000
```

5.4 NonDetects

The non detects in the Cq values and amplification efficiencies can be replaced by using the `NonDetects` function. This function is based on replacing the non detects within the technical replicates based on the mean or median values of the technical replicates of the sample. In the case of a sample, if all the technical replicates are present as non detects (NAs), then the non detects for that sample will not be replaced. It will replace the non detects in a sample, if atleast one of the technical replicate of a sample has expressed. The

`NonDetects` function works on the object of class `RTqPCRBatCh` and returns the resulting data too in the same format of class `RTqPCRBatCh`. The output of any of the functions, including `CqEfts`, `ReplaceAboveCutOff`, `ReplaceValue` can be implemented as an input for the `NonDetects` function. The default method of replacing the non detects is mean. Following is the example to replace the non detects on the basis of mean value within technical replicates and display of the first ten Cq and amplification efficiencies:

```
> Cqefts.nondetects <- NonDetects(res.LC)
> ## to express the overview of the resulting data
> Cqefts.nondetects

RTqPCRBatCh (storageMode: lockedEnvironment)
assayData: 96 features, 1 samples
  element names: efts, exprs
protocolData: none
phenoData: none
featureData
  featureNames: 1 *2* ... 96 (96 total)
  fvarLabels: Sample position Sample name ... sig.model (35 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

> ## to express the first ten Cq values
> exprs(Cqefts.nondetects)[1:10]

[1] 45.00    NA 24.87 29.37    NA 12.36  2.49    NA 20.69  2.49

> ## to express the first ten amplification efficiencies
> efts(Cqefts.nondetects)[1:10]

[1] 1.0024280          NA 0.6576489 1.0003550          NA 0.5977603 1.0000910
[8]          NA 0.5882534 1.0000910
```

Following is the example to replace the non detects on the basis of median value within technical replicates and display of the first ten Cq and amplification efficiencies:

```
> Cqefts.nondetects1 <- NonDetects(res.LC, Calc = "Median")
> ## to express the overview of the resulting data
> Cqefts.nondetects1
> ## to express the first ten Cq values
> exprs(Cqefts.nondetects1)[1:10]
> ## to express the first ten amplification efficiencies
> efts(Cqefts.nondetects1)[1:10]
```

Important Note: Before starting the further steps, we need to look on our final Cq values and amplification efficiencies, particularly of reference samples. The values of all the technical replicates of reference samples and positive calibrators of target and reference samples should not be NA, because reference samples are very important to compute the final fold concentration and relative expression ratio. If the final computed values of fold concentration and relative expression ratio has high number of NA, then the user needs to look the Cq values and amplification efficiencies of reference samples. The user can change the particular value by implementing the `ReplaceValue` function.

6 Combining the replicates, based on Cq values and amplification efficiencies

User can combine technical replicates of Cq values by default method (mean) and compute their standard deviations within replicates, through following code:

```
> ## Combine technical replicates of Cq values by default method (mean)
> Cqreps <- CombineTechReps(res.LC)
> ## to visualise the overview of the resulting data
> Cqreps
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
  element names:
protocolData: none
phenoData: none
featureData
  featureNames: 1 2 ... 30 (30 total)
  fvarLabels: ID Cq ... Combined sample and target type (5 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> ## to visualise the resulting data
> fData(Cqreps)[,c("ID", "Cq", "sd.Cq")]
```

	ID	Cq	sd.Cq
1	1	37.19	11.052
2	3	18.62	8.846
3	**7**	2.49	NA
4	9	20.93	0.332
5	15	18.44	0.245
6	21	18.50	NA


```

7  **22** 14.96 19.742
8      23 25.70      NA
9      27 12.60 16.405
10 **32**  7.94      NA
11      34  1.87      NA
12      36 13.43      NA
13 **37** 11.00      NA
14      42 20.09      NA
15      45 12.61      NA
16      51 17.96  8.224
17      54 11.33  6.527
18      58 17.67  0.594
19      60 19.55  3.437
20      63 19.30      NA
21      66 17.14  3.824
22      67 10.26  6.237
23 **74** 29.30      NA
24 **80** 10.48  0.396
25 **84** 26.82      NA
26      86 30.32  1.054
27      88 13.58      NA
28      91 12.28  3.745
29      94  6.91      NA
30 **95**  8.05      NA

```

```

> ## to visualise the resulting data along with the other information of replicates
> fData(Cqreps)

```

	ID	Cq	sd.Cq	Target name	Combined sample and target type
1	1	37.19	11.052	negativ Kontrolle	Target Negative
2	3	18.62	8.846	hk	Ref Unknown
3	**7**	2.49	NA	negativ Kontrolle	Target Negative
4	9	20.93	0.332	hk	Ref Unknown
5	15	18.44	0.245	hk	Ref Unknown
6	21	18.50	NA	hk	Ref Unknown
7	**22**	14.96	19.742	goi	Target PosCalibrator
8	23	25.70	NA	goi	Target Unknown
9	27	12.60	16.405	hk	Ref Unknown
10	**32**	7.94	NA	goi	Target Unknown
11	34	1.87	NA	goi	Target PosCalibrator
12	36	13.43	NA	hk	Ref Unknown

13	**37**	11.00	NA	goi	Target PosCalibrator
14	42	20.09	NA	hk	Ref Unknown
15	45	12.61	NA	hk	Ref Unknown
16	51	17.96	8.224	hk	Ref Unknown
17	54	11.33	6.527	hk	Ref Unknown
18	58	17.67	0.594	hk	Ref PosCalibrator
19	60	19.55	3.437	hk	Ref Unknown
20	63	19.30	NA	hk	Ref Unknown
21	66	17.14	3.824	hk	Ref Unknown
22	67	10.26	6.237	hk	Ref PosCalibrator
23	**74**	29.30	NA	goi	Target Unknown
24	**80**	10.48	0.396	goi	Ref Unknown
25	**84**	26.82	NA	hk	Ref Unknown
26	86	30.32	1.054	goi	Target Unknown
27	88	13.58	NA	goi	Target PosCalibrator
28	91	12.28	3.745	goi	Target Unknown
29	94	6.91	NA	goi	Target Unknown
30	**95**	8.05	NA	goi	Target Unknown

Combine technical replicates of amplification efficiencies and compute their standard deviation, by method of mean:

```
> ## Combine technical replicates of amplification efficiencies by default method (mean)
> Effsreps <- CombineTechReps(res.LC, cRepCq = FALSE)
> ## to visualise the overview of the resulting data
> Effsreps
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
  element names:
protocolData: none
phenoData: none
featureData
  featureNames: 1 2 ... 30 (30 total)
  fvarLabels: ID effs ... Combined sample and target type (7 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> ## to visualise the resulting data
> fData(Effsreps)[,c("ID", "effs", "sd.effs", "Cq", "sd.Cq")]
```

	ID	effs	sd.effs	Cq	sd.Cq
1	1	1.001	0.001466	37.19	11.052

2	3	0.628	0.042348	18.62	8.846
3	**7**	1.000	NA	2.49	NA
4	9	0.604	0.022614	20.93	0.332
5	15	0.760	0.293195	18.44	0.245
6	21	1.100	NA	18.50	NA
7	**22**	1.024	0.033736	14.96	19.742
8	23	2209.600	NA	25.70	NA
9	27	0.361	0.382880	12.60	16.405
10	**32**	0.281	NA	7.94	NA
11	34	0.879	NA	1.87	NA
12	36	0.663	NA	13.43	NA
13	**37**	0.359	NA	11.00	NA
14	42	1.096	NA	20.09	NA
15	45	0.511	NA	12.61	NA
16	51	0.915	0.367514	17.96	8.224
17	54	0.766	0.308475	11.33	6.527
18	58	0.822	0.381787	17.67	0.594
19	60	0.896	0.275894	19.55	3.437
20	63	1.097	NA	19.30	NA
21	66	0.675	0.065095	17.14	3.824
22	67	0.721	0.109550	10.26	6.237
23	**74**	0.978	NA	29.30	NA
24	**80**	0.988	0.000168	10.48	0.396
25	**84**	0.992	NA	26.82	NA
26	86	1.047	0.264689	30.32	1.054
27	88	1.132	NA	13.58	NA
28	91	1.047	0.118480	12.28	3.745
29	94	0.979	NA	6.91	NA
30	**95**	0.894	NA	8.05	NA

```
> ## to visualise the resulting data along with the other information of replicates
> fData(Effsreps)
```

	ID	effs	sd.effs	Cq	sd.Cq	Target name
1	1	1.001	0.001466	37.19	11.052	negativ Kontrolle
2	3	0.628	0.042348	18.62	8.846	hk
3	**7**	1.000	NA	2.49	NA	negativ Kontrolle
4	9	0.604	0.022614	20.93	0.332	hk
5	15	0.760	0.293195	18.44	0.245	hk
6	21	1.100	NA	18.50	NA	hk
7	**22**	1.024	0.033736	14.96	19.742	goi

8	23	2209.600	NA	25.70	NA	goi
9	27	0.361	0.382880	12.60	16.405	hk
10	**32**	0.281	NA	7.94	NA	goi
11	34	0.879	NA	1.87	NA	goi
12	36	0.663	NA	13.43	NA	hk
13	**37**	0.359	NA	11.00	NA	goi
14	42	1.096	NA	20.09	NA	hk
15	45	0.511	NA	12.61	NA	hk
16	51	0.915	0.367514	17.96	8.224	hk
17	54	0.766	0.308475	11.33	6.527	hk
18	58	0.822	0.381787	17.67	0.594	hk
19	60	0.896	0.275894	19.55	3.437	hk
20	63	1.097	NA	19.30	NA	hk
21	66	0.675	0.065095	17.14	3.824	hk
22	67	0.721	0.109550	10.26	6.237	hk
23	**74**	0.978	NA	29.30	NA	goi
24	**80**	0.988	0.000168	10.48	0.396	goi
25	**84**	0.992	NA	26.82	NA	hk
26	86	1.047	0.264689	30.32	1.054	goi
27	88	1.132	NA	13.58	NA	goi
28	91	1.047	0.118480	12.28	3.745	goi
29	94	0.979	NA	6.91	NA	goi
30	**95**	0.894	NA	8.05	NA	goi

Combined sample and target type

1	Target Negative
2	Ref Unknown
3	Target Negative
4	Ref Unknown
5	Ref Unknown
6	Ref Unknown
7	Target PosCalibrator
8	Target Unknown
9	Ref Unknown
10	Target Unknown
11	Target PosCalibrator
12	Ref Unknown
13	Target PosCalibrator
14	Ref Unknown
15	Ref Unknown
16	Ref Unknown
17	Ref Unknown

```

18             Ref PosCalibrator
19             Ref Unknown
20             Ref Unknown
21             Ref Unknown
22             Ref PosCalibrator
23             Target Unknown
24             Ref Unknown
25             Ref Unknown
26             Target Unknown
27             Target PosCalibrator
28             Target Unknown
29             Target Unknown
30             Target Unknown

```

Combine technical replicates based on Cq, by method of median:

```

> ## Combine technical replicates of Cq values by default method (mean)
> Cqreps1 <- CombineTechReps(res.LC, calc = "Median")
> ## to visualise the overview of the resulting data
> Cqreps1
> ## to visualise the resulting data
> fData(Cqreps1)[,c("ID", "Cq", "sd.Cq")]
> ## to visualise the resulting data along with the other information of replicates
> fData(Cqreps1)

```

Combine technical replicates based on amplification efficiencies, by method of median:

```

> ## Combine technical replicates of amplification efficiencies by default method (mean)
> Effsreps1 <- CombineTechReps(res.LC, calc = "Median", cRepCq = FALSE)
> ## to visualise the overview of the resulting data
> Effsreps1
> ## to visualise the resulting data
> fData(Effsreps1)[,c("ID", "effs", "sd.effs", "Cq", "sd.Cq")]
> ## to visualise the resulting data along with the other information of replicates
> fData(Effsreps1)

```

Combine technical replicates, based on Cq values by method of geometric mean:

```

> ## Combine technical replicates of Cq values by default method (mean)
> Cqreps2 <- CombineTechReps(res.LC, calc = "Geom")
> ## to visualise the overview of the resulting data
> Cqreps2
> ## to visualise the resulting data

```

```
> fData(Cqreps2)[,c("ID", "Cq", "sd.Cq")]
> ## to visualise the resulting data along with other information of replicates
> fData(Cqreps2)
```

Combine technical replicates, based on amplification efficiencies by method of geometric mean:

```
> ## Combine technical replicates of Cq values by Geometric mean method
> Effsreps2 <- CombineTechReps(res.LC, calc = "Geom", cRepCq = FALSE)
> ## to visualise the overview of the resulting data
> Effsreps2
> ## to visualise the resulting data
> fData(Effsreps2)[,c("ID", "effs", "sd.effs", "Cq", "sd.Cq")]
> ## to visualise the resulting data along with other information of replicates
> fData(Effsreps2)
```

7 Compute delta Cq, delta delta Cq values and fold concentration of the Target samples

This step is further divided into two parts: i) Computing delta Cq values of the combined technical replicates, based on the Cq values. ii) Computing the delta delta Cq values and log of fold concentration of the samples. Both the parts are illustrated in further subsections.

7.1 Compute delta Cq values

The delta Cq value can be computed by implementing the `DeltaCq` function on the combined Cq values of technical replicates. These combined Cq values can be produced from the `CombineTechReps` function. The object of class `RTqPCRBatch`, which is an output of `CombineTechReps` function, as the input and will populate an object of class `RTqPCRBatch` with new slots as an output. The following code can be implemented to calculate the delta Cq values:

```
> ## Combine technical replicates of Cq values by default method (mean)
> deltaCq <- DeltaCq(Cqreps, Ref = "hk")
> ## to visualise the overview of the resulting data
> deltaCq
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
  element names:
protocolData: none
phenoData: none
```

```

featureData
  featureNames: 8 10 ... 27 (11 total)
  fvarLabels: ID deltaCq ... Combined sample and target type (5 total)
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

> ## To visualise the results of delta Cq
> fData(deltaCq)

```

	ID	deltaCq	sd.deltaCq	Target	name	Combined sample and target type
8	23	8.03429	NA	goi		Target Unknown
10	**32**	-9.72571	NA	goi		Target Unknown
23	**74**	11.63429	NA	goi		Target Unknown
26	86	12.65429	6.07	goi		Target Unknown
28	91	-5.38571	7.06	goi		Target Unknown
29	94	-10.75571	NA	goi		Target Unknown
30	**95**	-9.61571	NA	goi		Target Unknown
7	**22**	-2.70571	20.63	goi		Target PosCalibrator
11	34	-15.79571	NA	goi		Target PosCalibrator
13	**37**	-6.66571	NA	goi		Target PosCalibrator
27	88	-4.08571	NA	goi		Target PosCalibrator

7.2 Compute delta delta Cq values and fold concentration

The `DeltaDeltaCqAll` function can be implemented to compute the delta delta Cq values and log of fold concentration of all the samples. The `DeltaDeltaCqAll` function will take an object of class `RTqPCRBatch` as an input and will populate an object of class `RTqPCRBatch` as an output. The output of `DeltaCq` function will act as the input for the function `DeltaDeltaCqAll`. The delta delta Cq values, standard deviation of delta delta Cq and fold concentration can be computed through following implementation:

```

> ## Compute delta delta Cq values and it's standard deviation and fold conentration
> deltadeltaCq <- DeltaDeltaCqAll(deltaCq)
> ## to visualise the overview of the resulting data
> deltadeltaCq
> ## to visualise the resulting data
> fData(deltadeltaCq)

```

8 Compute relative expression ratio

The function `NRQeffsAll` compute the relative expression ratio of Target by three methods of Roche, Pfaffl et al. and ddcq. The `NRQeffsAll` computes on the output of `CombineTechReps` function, where the technical replicates are combined on the basis of amplification efficiencies. Here is the example to implement the `NRQeffsAll` function:

```
> ## Compute delta delta Cq values and it's standard deviation and fold conentration
> exp.ratio <- NRQeffsAll(Effsreps, y= "hk")
> ## to visualise the overview of the resulting data
> exp.ratio
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
  element names:
protocolData: none
phenoData: none
featureData
  featureNames: 1 2 ... 7 (7 total)
  fvarLabels: ID Roche Method Pfaffl Method delta delta Cq Method
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:

> ## to visualise the resulting data
> fData(exp.ratio)
```

	ID	Roche	Method	Pfaffl	Method	delta	delta	Cq	Method
1	23		-199.269		-118.1844				10.64
2	**32**		8.714		-3.0624				-1.67
3	**74**		-0.713		0.4215				13.13
4	86		-2.757		-0.9171				13.84
5	91		-1.929		-0.0885				1.34
6	94		-1.218		-0.0731				-2.39
7	**95**		-0.463		-0.2580				-1.60

9 Auxiliary functions

This section is divided into five subsections, based on the functions:

9.1 read.Mx3005P

The `read.Mx3005P` function reads in the .txt Tab delimited files of raw fluorescence data of Mx3005P RTqPCR and populates an object of class `RTqPCR`. This function is used in the

read.RTqPCR function. An example to implement the read.Mx3005P function is as follows:

```
> path <- system.file("exData", package = "RTqPCR")
> Mx3005P.example <- file.path(path, "Mx3005P_Example.txt")
> rtData <- read.Mx3005P(Mx3005P.example)
> ##to express the overview of resulting data
> rtData
```

```
RTqPCRBatch (storageMode: lockedEnvironment)
assayData: 45 features, 34 samples
  element names: exprs
protocolData: none
phenoData
  sampleNames: 14 15 ... 78 (34 total)
  varLabels: Well Ramp/Plateau# ... mspSegment (5 total)
  varMetadata: labelDescription
featureData
  featureNames: 1 2 ... 45 (45 total)
  fvarLabels: Cycle# Temperature
  fvarMetadata: labelDescription
experimentData: use 'experimentData(object)'
Annotation:
```

```
> ## to express the phenoData of cycData
> head(pData(rtData))
```

	Well	Ramp/Plateau#	Dye	mspSegment
14	14	P	2 Campy	2
15	15	P	2 Campy	2
16	16	P	2 Campy	2
17	17	P	2 Campy	2
18	18	P	2 Campy	2
20	20	P	2 Campy	2

```
> ## to express the fData of cycData
> head(fData(rtData))
```

	Cycle#	Temperature
1	1	59.6
2	2	59.6
3	3	59.6
4	4	59.6
5	5	59.6
6	6	59.6

9.2 Read.LC480SampleInfo

The `Read.LC480SampleInfo` function reads in the .txt Tab delimited file of sample information data of LC480 light cycler and populate an object of class `AnnotatedDataFrame`. This function is a modified function from the `link[ReadqPCR]read.LC480SampleInfo`. One change is introduced in the previously designed function to read all columns, irrespective of the condition whether they are completely filled or empty. Reading all the columns is required to perform the computation of other functions of the work flow. Following is the example to implement the `Read.LC480SampleInfo` function:

```
> LC480.sampleInfo <- file.path(path, "LC480_example_SampleInfo.txt")
> LC480.saminfo <- Read.LC480SampleInfo(LC480.sampleInfo)
> ##To express the overview of sample information data
> LC480.saminfo
```

```
An object of class 'AnnotatedDataFrame'
 rowNames: 1 2 ... 96 (96 total)
 varLabels: Sample position Sample name ... Combined sample and target
            type (9 total)
 varMetadata: labelDescription
```

```
> ##To express the phenodata
> head(pData(LC480.saminfo))
```

	Sample position	Sample name	Replicate	of Filter combination	Target name
1	A1	Sample_1	A1	465-510	negativ Kontrolle
2	A2	Sample_2	A2	465-510	goi
3	A3	Sample_2	A3	465-510	hk
4	A4	Sample_3	A1	465-510	negativ Kontrolle
5	A5	Sample_4	A2	465-510	goi
6	A6	Sample_4	A3	465-510	hk

	Sample Pref color	Concentration	Efficiency	Combined sample and target type
1	\$00FF8000	NA	2	Target Negative
2	clRed	NA	2	Target Unknown
3	\$0030D700	NA	2	Ref Unknown
4	clFuchsia	NA	2	Target Negative
5	clGray	NA	2	Target Unknown
6	\$0012D7FA	NA	2	Ref Unknown

9.3 Read.Mx3005PSampleInfo

The `Read.Mx3005PSampleInfo` function reads in the .txt Tab delimited file of sample information data of Mx3005P RTqPCR and populate an object of class `AnnotatedDataFrame`.

This function is further used in `Read.RTqPCRSampleInfo` function. It's a supplementary function to `Read.RTqPCRSampleInfo` function. Following is the example to implement the `Read.Mx3005PSampleInfo` function:

```
> Mx3005P.sampleInfo <- file.path(path, "Mx3005P_example_SampleInfo.txt")
> Mx3005P.samInfo <- Read.Mx3005PSampleInfo(Mx3005P.sampleInfo)
> ##To express the overview of the sample information data
> Mx3005P.samInfo
```

An object of class 'AnnotatedDataFrame'

```
rowNames: 1 2 ... 34 (34 total)
varLabels: Well Replicate of Target name Combined sample and target
           type
varMetadata: labelDescription
```

```
> ##To express the phenodata
> head(pData(Mx3005P.samInfo))
```

	Well	Replicate	of Target name	Combined sample and target	type
1	14	A1	goi	Target	Unknown
2	15	A2	goi	Target	Unknown
3	16	A3	hk	Ref	Unknown
4	17	A1	goi	Target	Unknown
5	18	A2	goi	Target	Unknown
6	20	A3	hk	Ref	Unknown

9.4 DeltaDeltaCq

The `DeltaDeltaCq` function computes the delta delta Cq values and fold concentration of the sample. It is designed as a supplementary to `DeltaDeltaCqAll` function. The output of `DeltaCq` function acts as the input to this function. It acts on an object of class `RTqPCR-Batch` and returns the resulting data as an object of class `data.frame`. Although, it's an auxilliary function and can be directly implemented through `DeltaDeltaCqAll` function. Following is the example for implementing the `DeltaDeltaCq` function:

```
> x.deltadeltaCq <- DeltaDeltaCq(deltaCq)
> ## to express the resulting data
> x.deltadeltaCq
```

	ID	deltadeltaCq	Fold Conc.[log]	sd.deltadeltaCq
1	23	15.35	-10.64	NA
2	**32**	-2.41	1.67	NA

3	**74**	18.95	-13.14	NA
4	86	19.97	-13.84	21.5
5	91	1.93	-1.34	21.8
6	94	-3.44	2.38	NA
7	**95**	-2.30	1.59	NA

9.5 NRQeffs

The `NRQeffs` function computes the relative expression ratio and is designed as a supplementary function to the `NRQeffsAll` function. It takes the combined technical replicates, which are produced as an output of `CombineTechReps` function, where they are combined based on amplification efficiency. The `NRQeffs` function acts on an object of class `RTqPCRBatch` and returns the resulting data as an object of class `data.frame`. It's an auxiliary function to `NRQeffsAll`, which can be directly implemented through the `NRQeffs` function.

9.6 RTqPCR.dataframe

The `RTqPCR.dataframe` function is designed as a supplementary function to the graphical user interface (GUI) for package "RTqPCR". This GUI is based on the shinyApps. In shinyApps, the information are passed between functions as well as resulting ones are taken out in form of `data.frame`. So, to handle the functions based on S4 class `RTqPCRBatch`, the `RTqPCR.dataframe` function is designed, which can convert the object of class `data.frame` into object of class `RTqPCRBatch`. It also reveals a way to handle the S4 class object in shinyApp.

10 Graphical User Interface (GUI)

```
> #library(shiny)
> #RTqPCR.gui()
```